HAMMER

Body Screen TX


Art.-Nr.: 75000

HAMMER

LIEFERUMFANG
PACKAGE LIST

Body Screen TX

4*AA Batterien
4*AA Batteries

Gebrauchsanweisung
User Manual

Maßband
Tape

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01	Wichtige Sicherheitsrichtlinien	P. 01
02	Spezifikationen	P. 02
03	Lerne deine HAMMER Body Screen TX kennen	P. 03
04	Display	P. 04
05	App-Installation	P. 04
06	Verwendung der Body Screen TX	P. 05
07	Pflege & Wartung	P. 07
08	Garantiebedingungen	P. 07
09	Service-Kontakt	P. 08

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Table of content

01	Important Safety Guidelines	P. 10
02	Specifications	P. 11
03	Getting to Know Your Body Screen TX	P. 12
04	Display	P. 13
05	App Installation	P. 13
06	Using Your Body Screen TX	P. 14
07	Care & Maintenance	P. 16
08	Warranty	P. 16

Important Safety Guidelines

Read all instructions before using HAMMER BodyPedia Body Screen TX.

When using electrical appliance, basic safety precautions should always be followed, including the following:

- Read all instructions carefully before using this product. Retain this user manual for future reference.

WARNINGS

- Individuals with medical implant devices such as pacemakers, or essential support devices such as patient monitoring systems, must not use this equipment. Safe, low-level currents will flow through the body during the test, which may cause malfunctioning of the device or endanger lives.
- Bioelectrical Impedance Analysis (BIA) uses safe low level currents, which are not harmful to the body. However, we do not recommend pregnant women test.
- The scale surface may become slippery in humid places. Please watch your steps and don't slip.

CAUTION

- Stand upright for about 5 minutes before testing. Taking the test immediately after lying in bed or sitting for a long period of time might result in a slight change in the test results. This is because body water tends to move to the lower body as soon as the person stands or gets up.
- Do not eat before testing. In cases where the examinee has already eaten, the test should be put off for at least two hours after the meal. This is because food mass is included in the examinee's weight and thus, may result in measurement errors.
- Use the bathroom before testing. Waste is not included in the body's compositional elements, but the volume of urine and excrement is included in the weight measurement affecting accuracy of the test results.
- Do not exercise before testing. Strenuous exercise or sharp movements can cause temporary changes in body composition. Even light exercise can change your body composition temporarily.
- Take the test in the morning, if possible. Body water tends to gravitate towards the lower body throughout the day, affecting accuracy of the test results.
- Take the test at the same time and same condition in each day to reduce test variables, if possible, to achieve a coherent test record.
- Thoroughly wipe the palms and soles with wet tissue/towel before testing. Testing may be difficult if the examinee's palms and soles are too dry or if the examinee has too many calluses.
- Put the HAMMER Body Screen TX on a leveled, hard and non-vibrating surface. Put the scale on an uneven surface may cause the examinee to fall down. Test results may also be inaccurate.
- Avoid contact with the examinee during testing. Contact may lead to interference affecting test results.


Specifications

Bioelectrical Impedance Analysis (BIA) Measurement Items	Bioelectrical Impedance (z)	15 Impedance Measurements by Using 3 Different Frequencies (5kHz, 50kHz,250kHz)at Each of 5 Segment (Right Arm, Left Arm, Trunk, Right Leg, and Left Leg)
Electrode Method	Tetrapolar 8-Point Tactile Electrodes	
Measurement Method	Direct Segmental Multi-frequency Bioelectrical Impedance Analysis Method (DSM-BIA) Simultaneous Multi-frequency Impedance Measurement (SMFIM)	
Body Analyses	Weight, Body Fat Percentage, BMI, Bone Minerals, Total Body Water, Protein, Visceral Fat, Minerals, Muscle Mass, Body Fat Mass, Segmental Lean Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg), Waist-Hip-Ratio, Neck Width, Chest Width, Waist Width, Hip Width, Segmental Fat Analysis (Right Arm, Left Arm, Trunk, Right Leg, Left Leg), BodyPedia Score, Fat Control, Muscle Control, Impedance (Each segment and each frequency)	
Database	The firmware can save up to 10 results. The App can save up unlimited results.	
Connection	Bluetooth 5.0	
Battery Type	4 *AA Batteries	
Battery Life	Around 3 months	
Display Type	Hidden LED Screen	
Dimension	31.5x31.5x3.3cm	
Testing Time	About 20 seconds	
Weight	2.35kg	
Operation Environment	10~40°C	
Storage Environment	-10~60°C	
Testing Weight Range	5~200kg	
Testing Age Range	13 ~99 years	


Get To Know YOUR HAMMER BODY SCREEN TX

HAMMER Body Screen TX detects and analyzes body composition of different body parts by BIA (Bioelectrical Impedance Analysis). It uses 8-point tactile electrodes, multiple currents at low, medium and high frequencies, direct segmental measurement (directly measure each of the body's five cylinders left arm, right arm, torso, left leg, and right leg) and human body modeling to deliver quick, accurate and detailed results.

Product Diagram

DISPLAY

LED EFFECT	OFF	ON	Flash
	Not connected	Connected	Not connected
	Not measuring	Measuring	/
	Not connected	Connected	/
	Not connected	Connected	/
	Not measuring	Measuring	/

App Installation

Search "BodyPedia" in the Android Application Market, the IOS App Store or scan the following QR code to download and install the App

App downloaded via other channels may be unusable. Before pairing with your equipment, make sure your equipment meets the following requirements:

- iPhone 5S and above with iOS 10.0 and above system;
- Android 5.0 and above with Bluetooth 4.0

Scan QR Code
Download < BodyPedia App >

Connect to the smart scale via Bluetooth

Bind the scale to BodyPedia App


Open BodyPedia App and Bluetooth of the phone, follow the instruction of the App to register, bind the scale to the App and finish the first measurement.

Bluetooth connection


Step on the scale with one foot to activate the scale, and take your foot off the scale when the scale screen lights up. Open BodyPedia App, and the scale will connect to the App automatically if the scale has been bound to the App.

Using Your Body Screen TX

1. Pull the tab underneath the scale. Wait for the screen display 0.0kg/lb before weighing.

2. Step on the footplate barefoot, and weight measurement begins.

Stand Barefoot to Begin.

Put the scale on a leveled, non-vibrating surface.

3. Stand up straight for 3s, and weighing is completed when the figure stop flashing.




Figure stops flashing

4. Press and pull out the handles after weighing is completed.

5. The body test begins.

The examinee must maintain proper posture to have accurate test results.
 * The test will proceed when there is good electrical contact.

6. When the test is completed, connect the scale to your App to have the results.

Feet			
Legs			
Arms			
Hands			

Note: Please register in the BodyPedia App and follow the instructions to have your first test.

Care & Maintenance

Do not extend the handles beyond its limitation.

Do not place any objects on the footplate.

Do not apply excessive force on the equipment.

Do not allow any liquid substances to contact the equipment directly. Keep food and drinks away from the equipment. Substances getting inside the equipment can cause critical damage to the electronic components.

Remove the batteries when they become weak or when the body analyzer is not in use for extended periods. This will prevent damage to the device from battery leakage at such times

Use a lint-free cloth to gently wipe the external surface of the equipment about once every week

Be careful not to scratch the LED screen.

HAMMER Body Screen TX does not need regular maintenance. If some problems occur while operating the device, get in touch with our customer service team.

We do not take the responsibility about problems caused by any arbitrary repairs.

Warranty

We provide warranty for our devices according to the following conditions:

1) We will repair free of charge in accordance with the following conditions: Damage or defects to the device which are demonstrably due to a manufacturing defect, if they are reported to us immediately upon discovery and within 24 months after delivery to the end user. The warranty does not cover easily breakable parts, such as glass or plastic.

A warranty obligation is not triggered: by minor deviations from the target condition, which are insignificant for the value and usability of the device, by damage from chemical and electrochemical effects, from water and generally by damage due to force majeure.

2) The warranty service is provided in such a way that defective parts are repaired free of charge or replaced by faultless parts at our discretion. The costs for material and working time will be borne by us. Repairs at the place of installation cannot be demanded. The proof of purchase with the date of purchase and/or delivery must be presented. Replaced parts become our property.

3) The warranty claim expires if repairs or interventions are carried out by persons who are not authorized by us to do so or if our devices are equipped with supplementary or accessory parts which are not matched to our devices, furthermore if the device is damaged or destroyed by the influence of force majeure or by environmental influences, in the case of damage which has occurred due to improper handling - in particular non-observance of the operating instructions - or maintenance or if the device shows mechanical damage of any kind. The customer service can authorize you after telephone consultation to repair or exchange parts, which will be sent to you. In this case, the warranty claim will of course not expire.

4) Warranty services shall neither cause an extension of the warranty period nor shall they set a new warranty period in motion.

5) Further or other claims, in particular those for compensation for damage occurring outside the device, are excluded - insofar as liability is not mandatorily ordered by law.

6) Our warranty conditions, which include the prerequisites and scope of our warranty services, do not affect the contractual warranty obligations of the seller.

7) Parts subject to wear and tear, which are identified as such in the parts list, are not subject to the warranty conditions.

8) The warranty claim expires in case of improper use, especially in fitness studios, rehabilitation facilities and hotels.